

# Cutan® INSTANTGEL™

## ALCOHOL-BASED GEL HAND SANITISER

### PRODUCT DESCRIPTION

Highly effective broad spectrum alcohol based hand sanitising gel.

### WHERE TO USE

Ideal for use in all healthcare and food handling environments where high levels of hand hygiene are vital. Also, suitable for use in cafeterias, office areas, rest areas and any public facility to help reduce the spread of germs.

### BIOCIDAL EFFICACY

Bactericidal, Virucidal, Yeasticidal, Mycobactericidal.

### WHEN TO USE

**Healthcare environments** - use hand sanitiser as per the World Health Organisations '5 Moments of Hand Hygiene' recommendations.

**Food handling environments** - use hand sanitiser as recommended by HACCP risk assessment.

**General hygienic use** - use hand sanitiser after coughing, sneezing, contact with body fluids or touching potentially contaminated surfaces.

### HOW TO USE

Ensure hands are clean and dry before use. Use several times a day, as required.

**For general hygienic use:** apply 1-2 doses direct to clean hands and rub across all parts of dry and clean hands.

**For hygienic skin disinfection** (according to the EN1500 standard) - apply 3 doses (3ml) and spread across all parts of dry and clean hands for 30 seconds (60 seconds for viruses); re-apply as necessary to keep hands wet.

**For surgical disinfection:** (according to EN12791 standard) - apply 3 doses (3ml) and rub across all parts of dry and clean hands using standard hand rub procedure; re-apply as necessary to keep hands wet for 90 seconds.



	DESCRIPTION	BENEFITS
FEATURES	<b>Complete broad spectrum activity</b>	Highly effective formula kills 99.999% of many common germs, and is fully tested and proven to be bactericidal, virucidal, yeasticidal and mycobactericidal.
	<b>Full virucidal activity</b>	According to test method EN14476 the product is effective against all enveloped and non-enveloped viruses.
	<b>Instant gel action</b>	Rapidly and easily rubs across hands, providing the user with complete control of the product, preventing the dripping and splashing associated with using liquid alcohol sanitisers.
	<b>Pleasant-to-use</b>	The smooth gel format is designed for frequent use in-between hand washing to encourage hand hygiene compliance.
	<b>No water required</b>	Quickly dries on hands without the need for rinsing with water. Ideal for use outside the washroom at the point of need.
	<b>Skin Hypoallergenic</b>	This product is tested by an independent dermatologist using Human Repeat Insult Patch Tests (HRIPT) to verify and confirm this product does not irritate the skin and has very low allergenic potential.
INGREDIENTS	<b>Perfume-Free and Dye-Free</b>	Designed for people who are sensitive to perfumes and dyes and who prefer products which do not have these added.
	<b>80% alcohol content</b>	The formula with a high level of ethanol meets the requirements for effective broad-spectrum efficacy and is proven to be bactericidal, virucidal, yeasticidal, and mycobactericidal.
	<b>Contains Glycerin and Panthenol</b>	Skin moisturisers to help prevent drying and leave the skin feeling smooth after use.
	<b>Contains Hamamelis astringent extract</b>	Helps to prevent the skin from softening when sweating during glove use.
	<b>Silicone-free</b>	Reduces the risk of Silicone related compatibility issues, if in contact with surfaces prior to the application of paint or surface coatings.
	<b>Not made with ingredients of animal origin</b>	Helps customers with Halal and Kosher ingredient requirements.
COMPATIBILITY	<b>Suitable for food facilities</b>	The product has been independently evaluated and proven to not have the potential to taint food. Food handlers should use product only in accordance with usage instructions and product should not come directly into contact with food or items that contain food.
	<b>Proven glove compatibility</b>	Suitable for use with latex, vinyl and nitrile gloves.

# Cutan® INSTANTGEL™

## ALCOHOL-BASED GEL HAND SANITISER

### STATUTORY REGULATIONS

This product meets the current national regulations for Biocidal products.

### SAFETY DATA SHEETS

For Safety, Environmental, Handling, First Aid and Disposal information please refer to the Safety Data Sheet which can be downloaded from <https://www.scjp.com/en-gb/safety-data-sheets>

### HANDLING INSTRUCTIONS

Use biocides safely. Always read the label and product information before use. Please refer to the label or SDS for safety information.

### SHELF LIFE

This product has a shelf life of 36 months from the date of manufacture when stored unopened at room temperature.

### SUITABLE FOR USE IN FOOD HANDLING AREAS

When used as instructed, Cutan® INSTANTGEL™ does not influence the quality of foodstuffs. The product has been independently tested according to Sensory Analysis - Methodology - Triangle Test BS EN ISO 4120:2007 and has been proven to not have the potential to taint food, thereby supporting and meeting HACCP standards for impact on the smell and taste of foodstuffs.

### QUALITY ASSURANCE

SC Johnson Professional products are manufactured in facilities which follow appropriate quality assurance practices, including with respect to quality control processes for raw materials and finished goods, and comply with any regulatory requirements applicable to the manufacturing of such products.

### PRODUCT SAFETY ASSURANCE

#### Product Safety Assessment

The product has been assessed by an independent toxicologist and is declared safe for intended use on hands of adults and children above 3 years of age when used as directed, meeting all relevant regulatory requirements.

#### Skin Compatibility Test

Independent dermatological tests conducted to assess the skin compatibility of the product using voluntary test persons in semi occlusive, 48 consecutive hour application models, shows the product to be non-irritating in use.

#### Skin Hypoallergenic

The perfume-free and dye-free formulation did not demonstrate a potential for eliciting dermal irritation or sensitisation after conducting an independent 6 week Human Repeat Insult Patch Test (HRIPT).

Allergy sufferers should always refer to the ingredients list before using the product.

#### Glove Compatibility

In independent testing, the use of Cutan® INSTANTGEL™ hand sanitiser in combination with reputable latex, vinyl and nitrile gloves does not have any deleterious effect on the glove tensile properties using BS EN 455-2:2009 test methodology.

### PRODUCT TESTING

TEST TYPE	TEST STANDARD	TEST CONDITION	COMMENT
Bactericidal	EN1500	<i>In vivo</i>	Product passes at 30 seconds (3mls).
	EN12791	<i>In vivo</i>	Suitable for surgical hand disinfection – apply 3ml and keep hands wet for 90 seconds.
	EN1276	<i>In vitro</i>	Suitable for use in food, industrial, domestic and institutional areas. Possesses bactericidal activity in 30 seconds under clean conditions using 4 referenced strains.
	EN13727	<i>In vitro</i>	Suitable for use in the medical area. Possesses bactericidal activity in 30 seconds under clean conditions using 4 referenced strains.
Yeasticidal	EN1650	<i>In vitro</i>	Suitable for use in food, industrial, domestic and institutional areas. Possesses yeasticidal activity against <i>Candida albicans</i> at a contact time of 30 seconds.
	EN13624	<i>In vitro</i>	Suitable for use in the medical area. Possesses yeasticidal activity against <i>Candida albicans</i> at a contact time of 30 seconds.
Mycobactericidal	EN14348	<i>In vitro</i>	Suitable for use in the medical area. Possesses Mycobactericidal activity against <i>Mycobacterium terrae</i> and <i>Mycobacterium avium</i> in 30 seconds.
Virucidal	EN14476	<i>In vitro</i>	Deactivates Murine Norovirus (15 seconds), Adenovirus (30 seconds), Poliovirus (60 seconds) Virucidal according to EN14476 against all enveloped and non-enveloped viruses.

### INGREDIENTS

ALCOHOL DENAT., AQUA, GLYCERIN, PANTHENOL, ACRYLATES/C10-30 ALKYL ACRYLATE CROSSPOLYMER, HAMAMELIS VIRGINIANA (WITCH HAZEL) EXTRACT, HAMAMELIS VIRGINIANA (WITCH HAZEL) WATER, AMINOMETYL PROPANOL, TRISODIUM DICARBOXYMETHYL ALANINATE.

### PACK SIZES

STOCK CODE	SIZE	CASE QTY	DISPENSER
CAG50ML	50ml bottle	48	
CAG100ML	100ml Squeeze bottle	12	
CAG400P	400ml Pump bottle	12	
CAG392	1 Litre Cartridge	6	PROB01SA

SC Johnson Professional Limited  
Denby Hall Way  
Denby  
Derbyshire DE5 8JZ  
UK  
Tel: +44 (0) 1773 855100  
[www.scjp.com](http://www.scjp.com)

This information and all further technical advice is based upon our present knowledge and experience. However, it implies no liability or other legal responsibility on our part, including with regard to existing third party intellectual property rights, especially patent rights. In particular, no warranty, whether express or implied, or guarantee of product properties in the legal sense is intended or implied. We reserve the right to make any changes according to technological progress or further developments.

**scJohnson**  
PROFESSIONAL  
A family company®