



# Tork Couch Roll Advanced

150226

Think ahead.



## Description

Protect your patient couch with the disposable Tork Hygiene Roll Advanced. The perforated hygiene roll is ideal to maintain a high level of hygiene in a healthcare environment. Rolls can be used in Tork Couch Roll Dispenser.

- Better resistance and higher comfort due to 2-ply
- Perforation for easy tear off and consumption control
- Strong and easy to roll, for optimal protection
- Advanced
- Professional use

## Certifications



Tork  
Advanced

## Product Details

System	C1
--------	----

## Shipping Data

	Consumer Units (CON)	Transport unit (TRP)	Pallet (PAL)
EAN	7322541052403	7322541047782	7322541047799
Packaging Material	Plastic	Plastic	-
Pieces	2	18 (9 CON)	324 (18 TRP)
Height	500 mm	500 mm	1,161 mm
Length	120 mm	360 mm	1,200 mm
Width	120 mm	360 mm	1,000 mm
Gross Weight	1,023.97 g	9.3 kg	167.36 kg
Net Weight	980.1 g	8.82 kg	158.78 kg
Volume	7.2 dm <sup>3</sup>	64.8 dm <sup>3</sup>	1.39 m <sup>3</sup>
Layers Per Pallet	-	-	2
TRP Per Layer	-	-	9



Think ahead.

# Tork Couch Roll Advanced

150226

## Compatible Products



**Tork Disp Couch Roll Wall Stand**  
129183