



Manual Dispenser Swarfega®

PRODUCT DESCRIPTION

4L and 2L dispenser for sealed Swarfega® cartridges.

USAGE

The dispenser is for exclusive use with the Swarfega® range of 4L and 2L cartridges.

METHOD OF USE

Push the large button fully to dispense a measured dose of product appropriate for the application required. Empty cartridges are easily replaced in seconds.



4L Dispenser



2L Dispenser

	DESCRIPTION	BENEFITS
FEATURES	Guaranteed for Life	Swarfega® 4L and 2L Dispenser is guaranteed for lifetime. Guarantee applies only in the European Union, United Kingdom, and Ireland and to dispensers exclusively using SC Johnson Professional® products. Other terms and conditions apply: https://www.scjp.com/en-gb/guarantee
	Buttons protected by BioCote® ¹ technology	BioCote® silver-ion technology effectively reduces bacteria, mould and fungi on the surface of BioCote® treated SC Johnson Professional® dispenser buttons within as little as 2 hours and achieves up to 99.99% reduction over a 24 hour period*.
	Large easy-push buttons	Push the large button fully to dispense a measured dose of product appropriate for the application required.
	Measured dose	Delivers the right amount for effective product application, controlling over-usage and waste
	Hygienic, airless cartridge	The ultrasonically sealed cartridge and pump mechanisms allow no ingress of air into the cartridge during use to help prevent product contamination and ensure 'fresh' product is dispensed every time.
	Minimal waste	The air-tight cartridges are designed to collapse during use to minimise product waste.
HANDLING	Secure lockable mechanism	Unique mechanism helps prevent unauthorised access, pilferage and product contamination.
	Sight glass	Indicates the presence of a cartridge.
	Time saving	Cartridges can be replaced in seconds, minimising maintenance time and associated costs.

STOCK CODE	DISPENSER	QTY PER BOX	DIMENSIONS	APPLICATIONS PER CARTRIDGE	PRODUCT COMPATIBILITY/ INSERT	DOSE SIZE
SWA4000D	Swarfega® 4L Dispenser	6	348H x 177W x 160Dmm	888 doses	4L cartridges - Blue Insert Swarfega® Orange-SORC4LTR Swarfega® EXTREME-SHD4LTR	4.5ml
SWA2000D	Swarfega® 2L Dispenser	8	289H x 162W x 138Dmm	444 doses	2L Cartridges - Blue Insert Swarfega® LIGHT - SWC2LTSC Swarfega® Orange - SOR2LT	4.5ml

¹ BioCote® is a registered trademark of Biocote Limited *BioCote® maintains its antimicrobial protection for the lifetime of the dispenser button and will not wear off, wash off or leach out. BioCote® technology is not a substitute for good hand hygiene practices.

SC Johnson Professional Limited
Denby Hall Way
Denby
Derbyshire DE5 8JZ
UK
Tel: +44 (0) 1773 855100
www.scjp.com

This information and all further technical advice is based upon our present knowledge and experience. However, it implies no liability or other legal responsibility on our part, including with regard to existing third party intellectual property rights, especially patent rights. In particular, no warranty, whether express or implied, or guarantee of product properties in the legal sense is intended or implied. We reserve the right to make any changes according to technological progress or further developments.

© 2025 SC Johnson Professional Limited